



REF. DOCTOR : doctor@example.com

PATIENT : Recovery Test | ID: #100100

AGE/SEX : 30 Years
DRAWN : 28/03/2026 17:48:15
RECEIVED : 28/03/2026 17:48:15
REPORTED : 28/03/2026 17:48:15

Test Report Status	Final	Results	Reference Interval	Units
AI-BASED TUBERCULOSIS DIAGNOSTIC ANALYSIS (TBFERON-AI)				
SPECIMEN SOURCE METHOD : DIGITAL CHEST X-RAY + CLINICAL DATA		MULTIMODAL INPUT		
CNN X-RAY PROBABILITY METHOD : DEEP LEARNING - MOBILENETV2		52.9%	< 40.0%	%
CLINICAL PARAMETER SCORE METHOD : LOGISTIC REGRESSION ML MODEL		5.0%	< 40.0%	%
FINAL FUSED TB PROBABILITY METHOD : WEIGHTED MULTIMODAL FUSION (60% IMAGE + 40% CLINICAL)		33.7% Low	< 40.0%	%
COUGH DURATION METHOD : CLINICAL PARAMETER		1 Weeks	> 2 Weeks indicates risk	Weeks
FEVER / WEIGHT LOSS / NIGHT SWEATS METHOD : CLINICAL PARAMETER		No / No / No		
SPUTUM SMEAR MICROSCOPY METHOD : LAB DATA / ACID-FAST BACILLI		Not Done	NEGATIVE	
GENEXPERT MTB/RIF METHOD : LAB DATA / NUCLEIC ACID AMPLIFICATION		Not Done	NEGATIVE	
INTERPRETATION METHOD : AI MULTIMODAL FUSION ANALYSIS		NEGATIVE	NEGATIVE	

LUNG REGION CLASSIFICATION (MedSAM AI Segmentation)		
Anatomical Region	Severity	Status
Upper Right Lobe	HIGH	Critically Affected
Upper Left Lobe	LOW	Mildly Affected
Mid Right Lobe	HIGH	Critically Affected
Mid Left Lobe	LOW	Mildly Affected
Lower Lobes	HIGH	Critically Affected

Interpretation(s)

The AI-based TB detection system uses a multimodal approach combining Deep Learning analysis of the Chest X-Ray image with clinical parameter scoring (cough duration, fever, weight loss, night sweats) to arrive at a fused probability score.

AI Recommendation: TB unlikely. Consider other respiratory conditions if symptoms persist.

DISCLAIMER: This report is generated by an AI system for clinical decision support only. It is not intended to replace professional medical diagnosis. A qualified physician should evaluate these results in the context of the patient's full clinical picture. Confirmatory tests (sputum culture, GeneXpert, TST) are recommended for definitive TB diagnosis.

Authorized Signatory
TB-Vision AI Diagnostic System

Reviewing Physician
doctor@example.com



