

Patient Photo

REF. DOCTOR : doctor@example.com

PATIENT : pravin | ID: #958517

AGE/SEX : 49 Years
DRAWN : 01/04/2026 12:59:29
RECEIVED : 01/04/2026 12:59:29
REPORTED : 01/04/2026 12:59:29

Test Report Status	Final	Results	Reference Interval	Units
AI-BASED TUBERCULOSIS DIAGNOSTIC ANALYSIS (TBFERON-AI)				
SPECIMEN SOURCE METHOD : DIGITAL CHEST X-RAY + CLINICAL DATA		MULTIMODAL INPUT		
CNN X-RAY PROBABILITY METHOD : DEEP LEARNING - MOBILENETV2		78.1%	< 40.0%	%
CLINICAL PARAMETER SCORE METHOD : LOGISTIC REGRESSION ML MODEL		99.0%	< 40.0%	%
FINAL FUSED TB PROBABILITY METHOD : WEIGHTED MULTIMODAL FUSION (60% IMAGE + 40% CLINICAL)		86.4% High	< 40.0%	%
COUGH DURATION METHOD : CLINICAL PARAMETER		9 Weeks	> 2 Weeks indicates risk	Weeks
FEVER / WEIGHT LOSS / NIGHT SWEATS METHOD : CLINICAL PARAMETER		Yes / Yes / No		
SPUTUM SMEAR MICROSCOPY METHOD : LAB DATA / ACID-FAST BACILLI		Positive	NEGATIVE	
GENEXPERT MTB/RIF METHOD : LAB DATA / NUCLEIC ACID AMPLIFICATION		Negative	NEGATIVE	
INTERPRETATION METHOD : AI MULTIMODAL FUSION ANALYSIS		POSITIVE	NEGATIVE	

LUNG REGION CLASSIFICATION (MedSAM AI Segmentation)		
Anatomical Region	Severity	Status
Upper Right Lobe	MEDIUM	Moderately Affected
Upper Left Lobe	HIGH	Critically Affected
Mid Right Lobe	LOW	Mildly Affected
Mid Left Lobe	MEDIUM	Moderately Affected
Lower Lobes	MEDIUM	Moderately Affected

Interpretation(s)

The AI-based TB detection system uses a multimodal approach combining Deep Learning analysis of the Chest X-Ray image with clinical parameter scoring (cough duration, fever, weight loss, night sweats) to arrive at a fused probability score.

AI Recommendation: Immediate medical evaluation required. Isolate patient and schedule sputum test.

DISCLAIMER: This report is generated by an AI system for clinical decision support only. It is not intended to replace professional medical diagnosis. A qualified physician should evaluate these results in the context of the patient's full clinical picture. Confirmatory tests (sputum culture, GeneXpert, TST) are recommended for definitive TB diagnosis.

Authorized Signatory
TB-Vision AI Diagnostic System

Reviewing Physician
doctor@example.com



